

## CAB20676

---

### Product Information

<b>Product SKU:</b>	CAB20676	<b>Gene ID:</b>	4331372	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Oryza sativa

---

### Additional Information

<b>Observed MW:</b>	100kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	20kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Ehd4 is a novel Oryza-genus-specific regulator of Ehd1, and it plays an essential role in photoperiodic control of flowering time in rice.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** Ehd4; OsC3H20

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of Oryza sativa Ehd4.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.