

CAB20677

Product Information

Product SKU:	CAB20677	Gene ID:	4331345	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	92kDa	Conjugate:	Unconjugated
Calculated MW:	92kDa	Isotype:	IgG

Immunogen Information

Background: the mutation of LF1 induces ectopic expression of OSH1, which results in the initiation of lateral meristems to generate lateral florets in the axil of the sterile lemma.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: hb1; HOX10; OsHB1; OsHox10; OsJ_008792; LF1

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa LF1.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.