

CAB20686

Product Information

Product SKU:	CAB20686	Gene ID:	60343	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	25kDa	Conjugate:	Unconjugated
Calculated MW:	25kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a cytokine-like protein. The expression of this gene may be regulated by peroxisome proliferator-activated receptor gamma, and the encoded protein may be involved in the regulation of glucose and lipid metabolism. Alternative splicing results in multiple transcript variants.

Recommended Dilution: WB, 1:500 - 1:2000 IHC-P, 1:50 - 1:200

Synonyms: DLD; 2.19; XAP-7; DXS560S; FAM3A

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 34-230 of human FAM3A (NP_068578.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.