

CAB20692

Product Information

Product SKU: CAB20692	Gene ID: 4338011	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Oryza sativa

Additional Information

Observed MW: 53kDa	Conjugate: Unconjugated
Calculated MW: 49kDa	Isotype: IgG

Immunogen Information

Background: GW5 acts in the brassinosteroid signalling pathway to regulate grain width and weight in rice

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: GW5

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa GW5.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.