

CAB20703

Product Information

Product SKU: CAB20703	Gene ID: 4331799	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Oryza sativa

Additional Information

Observed MW: 37kDa	Conjugate: Unconjugated
Calculated MW: 33kDa	Isotype: IgG

Immunogen Information

Background: Involved in auxin transport. Positive regulator of the auxin signaling pathway involved in gravitropism, plant growth and grain development.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: BG1; OsJ_09620

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa BG1.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.