

## CAB20704

---

### Product Information

<b>Product SKU:</b>	CAB20704	<b>Gene ID:</b>	4346677	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Oryza sativa

---

### Additional Information

<b>Observed MW:</b>	70kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	45kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** The Histone Methyltransferase SDG724 Mediates H3K36me2/3 Deposition at MADS50 and RFT1 and Promotes Flowering in Rice.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** SDG724; lvp1; OsSET34

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of Oryza sativa SDG724.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.