

CAB20705

Product Information

Product SKU:	CAB20705	Gene ID:	4337864	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	70kDa	Conjugate:	Unconjugated
Calculated MW:	79kDa	Isotype:	IgG

Immunogen Information

Background: Chalk5 encodes a vacuolar H⁺-translocating pyrophosphatase influencing grain chalkiness in rice.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: Chalk5

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa Chalk5.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.