

CAB20713

Product Information

Product SKU: CAB20713	Gene ID: 4336878	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Oryza sativa

Additional Information

Observed MW: 53kDa	Conjugate: Unconjugated
Calculated MW: 53kDa	Isotype: IgG

Immunogen Information

Background: Involved in chilling tolerance. Interacts with the G-protein alpha subunit GPA1 to activate the calcium channel for sensing low temperature and to accelerate G-protein GTPase activity.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: COLD1; OsJ_16035; Cold1

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa Cold1.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.