

## CAB20715

---

### Product Information

<b>Product SKU:</b> CAB20715	<b>Gene ID:</b> 4330957	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Oryza sativa

---

### Additional Information

<b>Observed MW:</b> 40kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 40kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** OsMKK4 is a factor for grain size, can influence brassinosteroid (BR) responses and the expression of BR-related genes.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** OJ1316\_E06.24; OsMKK4

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of Oryza sativa OsMKK4.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.