

CAB20728

Product Information

| | | | | | |
|---------------------|----------|----------------------|--------|--------------------|-------------|
| Product SKU: | CAB20728 | Gene ID: | 79792 | Size: | 20uL, 100uL |
| Clone No: | ARC50993 | Host Species: | Rabbit | Reactivity: | Human |

Additional Information

| | | | |
|-----------------------|-------------|-------------------|--------------|
| Observed MW: | 53kDa/22kDa | Conjugate: | Unconjugated |
| Calculated MW: | 53kDa | Isotype: | IgG |

Immunogen Information

Background: Gasdermin D is a member of the gasdermin family. Members of this family appear to play a role in regulation of epithelial proliferation. Gasdermin D has been suggested to act as a tumor suppressor. Alternatively spliced transcript variants have been described.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: DF5L; DFNA5L; FKSG10; GSDMDC1; GSDMD

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 346-484 of human GSDMD (NP_079012.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH 7.3.