

CAB20733

Product Information

Product SKU:	CAB20733	Gene ID:	829389	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Arabidopsis thaliana

Additional Information

Observed MW:	46kDa	Conjugate:	Unconjugated
Calculated MW:	46kDa	Isotype:	IgG

Immunogen Information

Background: Mutant has elevated levels of free IAA in dominant mutant allele; Flavin Monooxygenase-Like Enzyme; Auxin Biosynthesis

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: L23H3.20; L23H3_20; YUC; YUCCA; YUCCA 1; YUCCA1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 71-221 of arabidopsis thaliana YUCCA1 (NP_194980.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.