

CAB20738

Product Information

Product SKU:	CAB20738	Gene ID:	114897	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	32kDa	Conjugate:	Unconjugated
Calculated MW:	32kDa	Isotype:	IgG

Immunogen Information

Background: Enables collagen binding activity. Involved in several processes, including negative regulation of platelet aggregation; positive regulation of aldosterone secretion; and positive regulation of cytosolic calcium ion concentration. Located in extracellular space. Is integral component of plasma membrane.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: GIP; CTRP1; ZSIG37; C1QTNF1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 142-281 of human C1QTNF1 (NP_112230.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.