

CAB20749

Product Information

Product SKU: CAB20749	Gene ID: 4342431	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Oryza sativa

Additional Information

Observed MW: 138kDa	Conjugate: Unconjugated
Calculated MW: 138kDa	Isotype: IgG

Immunogen Information

Background: Central factor in ethylene signaling pathways that control development, senescence and grain size. Acts as positive component of the ethylene-signaling pathway.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: EIN2; MHZ7; OsEIN2; P0455F03.4

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa OsEIN2

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.