

CAB20757

Product Information

Product SKU: CAB20757	Gene ID: 92255	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Mouse

Additional Information

Observed MW: 81kDa	Conjugate: Unconjugated
Calculated MW: 81kDa	Isotype: IgG

Immunogen Information

Background: Involved in adrenergic receptor signaling pathway. Located in membrane.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: DENBA; LMBRD2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 210-390 of human LMBRD2 (NP_001007528.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.