

## CAB20758

---

### Product Information

<b>Product SKU:</b> CAB20758	<b>Gene ID:</b> 8287	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b> 291kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 291kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** This gene is a member of the peptidase C19 family. It encodes a protein that is similar to ubiquitin-specific proteases, which cleave the ubiquitin moiety from ubiquitin-fused precursors and ubiquitinated proteins.

**Recommended Dilution:** WB,1:100 - 1:500 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

**Synonyms:** DFFRY; SPGFY2; USP9X/USP9Y

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 2456-2555 of human USP9X/USP9Y (NP\_004645.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.