

CAB20766

Product Information

Product SKU:	CAB20766	Gene ID:	7077	Size:	20uL, 100uL
Clone No:	ARC51187	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat,African green monkey

Additional Information

Observed MW:	22kDa	Conjugate:	Unconjugated
Calculated MW:	24kDa	Isotype:	IgG

Immunogen Information

Background: This gene is a member of the TIMP gene family. The proteins encoded by this gene family are natural inhibitors of the matrix metalloproteinases, a group of peptidases involved in degradation of the extracellular matrix. In addition to an inhibitory role against metalloproteinases, the encoded protein has a unique role among TIMP family members in its ability to directly suppress the proliferation of endothelial cells. As a result, the encoded protein may be critical to the maintenance of tissue homeostasis by suppressing the proliferation of quiescent tissues in response to angiogenic factors, and by inhibiting protease activity in tissues undergoing remodelling of the extracellular matrix.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: DDC8; CSC-21K; TIMP2

Purification Method: Affinity purification

Immunogen: Recombinant protein of human TIMP2.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.