

## CAB20768

---

### Product Information

<b>Product SKU:</b>	CAB20768	<b>Gene ID:</b>	8404	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	140kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	75kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Predicted to enable calcium ion binding activity; collagen binding activity; and extracellular matrix binding activity. Predicted to be involved in anatomical structure development and regulation of synapse organization. Located in extracellular space.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** SC1; MAST9; PIG33; MAST 9; SPARCL1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 17-300 of human SPARCL1 (NP\_004675.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.