

CAB20776

Product Information

Product SKU:	CAB20776	Gene ID:	10084	Size:	20uL, 100uL
Clone No:	ARC2697	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	38kDa	Conjugate:	Unconjugated
Calculated MW:	30kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a nuclear polyglutamine-binding protein that is involved with transcription activation. The encoded protein contains a WW domain. Mutations in this gene have been found in patients with Renpenning syndrome 1 and other syndromes with X-linked cognitive disability. Multiple alternatively spliced transcript variants that encode different protein isoforms have been described for this gene.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: SHS; MRX2; MRX55; MRXS3; MRXS8; NPW38; RENS1; PQBP1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-265 of human PQBP1 (NP_001027553.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.