

## CAB20802

---

### Product Information

<b>Product SKU:</b>	CAB20802	<b>Gene ID:</b>	105	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	81kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	81kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a member of the double-stranded RNA adenosine deaminase family of RNA-editing enzymes and may play a regulatory role in RNA editing.

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** RED2; ADAR3; ADARB2

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 600-739 of human ADARB2 (NP\_061172.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.