

CAB20804

Product Information

Product SKU:	CAB20804	Gene ID:	137964	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	50kDa	Conjugate:	Unconjugated
Calculated MW:	52kDa	Isotype:	IgG

Immunogen Information

Background: Lysophosphatidic acid acyltransferases (EC 2.3.1.51) catalyze the conversion of lysophosphatidic acid (LPA) to phosphatidic acid (PA). LPA and PA are involved in signal transduction and lipid biosynthesis.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: AGPAT6; LPAATZ; TSARG7; 1-AGPAT 6; LPAAT-zeta; GPAT4

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 205-456 of human GPAT4 (NP_848934.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.