

CAB20809

Product Information

| | | | | | |
|---------------------|----------|----------------------|--------|--------------------|-------------|
| Product SKU: | CAB20809 | Gene ID: | 91607 | Size: | 20uL, 100uL |
| Clone No: | ARC51651 | Host Species: | Rabbit | Reactivity: | Human |

Additional Information

| | | | |
|-----------------------|--------|-------------------|--------------|
| Observed MW: | 101kDa | Conjugate: | Unconjugated |
| Calculated MW: | 103kDa | Isotype: | IgG |

Immunogen Information

Background: Enables tRNA binding activity. Involved in several processes, including defense response to virus; negative regulation of G1/S transition of mitotic cell cycle; and replication fork arrest. Located in cytosol; nucleoplasm; and site of DNA damage.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: SLFN8/9; SLFN11

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 132-328 of human SLFN11 (NP_689483.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.