

CAB20840

Product Information

Product SKU:	CAB20840	Gene ID:	7421	Size:	20uL, 100uL
Clone No:	ARC51248	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	48kDa	Conjugate:	Unconjugated
Calculated MW:	48kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes vitamin D3 receptor, which is a member of the nuclear hormone receptor superfamily of ligand-inducible transcription factors. This receptor also functions as a receptor for the secondary bile acid, lithocholic acid. Downstream targets of vitamin D3 receptor are principally involved in mineral metabolism, though this receptor regulates a variety of other metabolic pathways, such as those involved in immune response and cancer. Mutations in this gene are associated with type II vitamin D-resistant rickets. A single nucleotide polymorphism in the initiation codon results in an alternate translation start site three codons downstream. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. A recent study provided evidence for translational readthrough in this gene, and expression of an additional C-terminally extended isoform via the use of an alternative in-frame translation termination codon.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: NR111; PPP1R163; VDR

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human VDR (NP_000367.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.