

CAB20841

Product Information

Product SKU:	CAB20841	Gene ID:	311	Size:	20uL, 100uL
Clone No:	ARC51481	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	56kDa	Conjugate:	Unconjugated
Calculated MW:	54kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the annexin family, a group of calcium-dependent phospholipid-binding proteins. Annexins have unique N-terminal domains and conserved C-terminal domains, which contain calcium-dependent phospholipid-binding sites. The encoded protein is a 56-kD antigen recognized by sera from patients with various autoimmune diseases. Several transcript variants encoding two different isoforms have been identified.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: ALS23; ANX11; CAP50; CAP-50; IBMWMA; Annexin A11

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 276-505 of human Annexin A11 (NP_001148.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.