

CAB20861

Product Information

Product SKU:	CAB20861	Gene ID:	546644	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	20-40kDa	Conjugate:	Unconjugated
Calculated MW:	14kDa	Isotype:	IgG

Immunogen Information

Background: Located in external side of plasma membrane.

Recommended Dilution: WB, 1:500 - 1:1000 IF/ICC, 1:50 - 1:200

Synonyms: Gr1; Gr-1; Ly-6G

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 27-105 of mouse Ly-6G (NP_001297367.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.