

CAB20866

Product Information

Product SKU:	CAB20866	Gene ID:	3948	Size:	20uL, 100uL
Clone No:	ARC2796	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	36kDa	Conjugate:	Unconjugated
Calculated MW:	36kDa	Isotype:	IgG

Immunogen Information

Background: Lactate dehydrogenase C catalyzes the conversion of L-lactate and NAD to pyruvate and NADH in the final step of anaerobic glycolysis. LDHC is testis-specific and belongs to the lactate dehydrogenase family. Two transcript variants have been detected which differ in the 5' untranslated region.

Recommended Dilution: WB, 1:100 - 1:500

Synonyms: CT32; LDH3; LDHX; Lactate Dehydrogenase C

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Lactate Dehydrogenase C (P07864).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.