

## CAB20870

---

### Product Information

<b>Product SKU:</b>	CAB20870	<b>Gene ID:</b>	15896	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2802	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	45kDa/55kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	31kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Predicted to enable integrin binding activity. Predicted to be involved in cell adhesion. Located in microvillus; plasma membrane; and uropod. Is integral component of plasma membrane. Is expressed in central nervous system; extraembryonic component; female reproductive system; retina; and vascular element. Orthologous to human ICAM2 (intercellular adhesion molecule 2).

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** CD102; Ly-60; Icam-2; ICAM2 / CD102

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of mouse ICAM2 / CD102 (P35330).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.