

CAB20891

Product Information

Product SKU:	CAB20891	Gene ID:	9322	Size:	20uL, 100uL
Clone No:	ARC2844	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	76kDa	Conjugate:	Unconjugated
Calculated MW:	68kDa	Isotype:	IgG

Immunogen Information

Background: Enables identical protein binding activity. Predicted to be involved in actin cytoskeleton organization; endocytosis; and signal transduction. Located in nucleoplasm. Biomarker of Huntington's disease.

Recommended Dilution: WB,1:100 - 1:500

Synonyms: STP; CIP4; HSTP; STOT; TRIP-10; Cip4

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cip4 (Q15642).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.