

## CAB20904

---

### Product Information

<b>Product SKU:</b>	CAB20904	<b>Gene ID:</b>	9770	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2868	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	38kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	38kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a protein that contains a Ras association domain. Similar to its cattle and sheep counterparts, this gene is located near the prion gene. Two alternatively spliced transcripts encoding the same isoform have been reported.

**Recommended Dilution:** WB,1:100 - 1:500

**Synonyms:** CENP-34; RASFADIN; RASSF2

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 227-326 of human RASSF2 (P50749).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.