

CAB20919

Product Information

Product SKU:	CAB20919	Gene ID:	5303	Size:	20uL, 100uL
Clone No:	ARC2893	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	14kDa	Conjugate:	Unconjugated
Calculated MW:	14kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the parvulin subfamily of the peptidyl-prolyl cis/trans isomerase protein family. The encoded protein catalyzes the isomerization of peptidylprolyl bonds, and may play a role in the cell cycle, chromatin remodeling, and/or ribosome biogenesis. The encoded protein may play an additional role in the mitochondria.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: EPVH; PAR14; PAR17; hEPVH; hPar14; hPar17; PIN4

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 52-131 of human PIN4 (Q9Y237).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.