

## CAB20923

---

### Product Information

<b>Product SKU:</b>	CAB20923	<b>Gene ID:</b>	56521	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2898	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse

---

### Additional Information

<b>Observed MW:</b>	23kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	23kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a member of a subclass of the HSP40/DnaJ protein family. Members of this family of proteins are associated with complex assembly, protein folding, and export. Two transcript variants encoding distinct isoforms have been identified for this gene.

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** JDP1; HPANBH4; DNAJC12

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 99-198 of human DNAJC12 (Q9UKB3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.