

## CAB20925

---

### Product Information

<b>Product SKU:</b>	CAB20925	<b>Gene ID:</b>	11336	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2900	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	85kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	86kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** The protein encoded by this gene is a component of the exocyst complex, a multiple protein complex essential for targeting exocytic vesicles to specific docking sites on the plasma membrane. Though best characterized in yeast, the component proteins and functions of exocyst complex have been demonstrated to be highly conserved in higher eukaryotes. At least eight components of the exocyst complex, including this protein, are found to interact with the actin cytoskeletal remodeling and vesicle transport machinery. The complex is also essential for the biogenesis of epithelial cell surface polarity.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** SEC6; Sec6p; SEC6L1; EXOC3

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 650-745 of human EXOC3 (O60645).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.