

CAB20929

Product Information

Product SKU:	CAB20929	Gene ID:	9672	Size:	20uL, 100uL
Clone No:	ARC2907	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	70kDa	Conjugate:	Unconjugated
Calculated MW:	45kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene belongs to the syndecan proteoglycan family. It may play a role in the organization of cell shape by affecting the actin cytoskeleton, possibly by transferring signals from the cell surface in a sugar-dependent mechanism. Allelic variants of this gene have been associated with obesity.

Recommended Dilution: WB,1:100 - 1:500

Synonyms: SDCN; SYND3; Syndecan 3

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Syndecan 3 (O75056).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.