

CAB20935

Product Information

Product SKU:	CAB20935	Gene ID:	285237	Size:	20uL, 100uL
Clone No:	ARC2917	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	38kDa	Conjugate:	Unconjugated
Calculated MW:	38kDa	Isotype:	IgG

Immunogen Information

Background: Involved in positive regulation of apoptotic process. Located in nucleus.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: -

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 51-200 of human C3orf38 (Q5JPI3).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.