

CAB20937

Product Information

Product SKU:	CAB20937	Gene ID:	80179	Size:	20uL, 100uL
Clone No:	ARC2919	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat,Monkey

Additional Information

Observed MW:	109kDa	Conjugate:	Unconjugated
Calculated MW:	109kDa	Isotype:	IgG

Immunogen Information

Background: Enables actin binding activity. Involved in regulation of cytokinesis and regulation of mitochondrial fission. Acts upstream of or within mitochondrion migration along actin filament. Located in cytosol and mitochondrial outer membrane.

Recommended Dilution: WB,1:100 - 1:500

Synonyms: MYOHD1; MYO19

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 450-650 of human MYO19 (Q96H55).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.