

CAB20953

Product Information

Product SKU:	CAB20953	Gene ID:	8799	Size:	20uL, 100uL
Clone No:	ARC2940	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	28kDa	Conjugate:	Unconjugated
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene facilitates peroxisomal proliferation and interacts with PEX19. The encoded protein is found in the peroxisomal membrane. Several transcript variants, some protein-coding and some not protein-coding, have been found for this gene.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: PEX14B; PEX11beta; PEX11-BETA; PEX11B

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human PEX11B (O96011).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.