

## CAB20955

---

### Product Information

<b>Product SKU:</b>	CAB20955	<b>Gene ID:</b>	57122	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2942	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Rat

---

### Additional Information

<b>Observed MW:</b>	106kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	106kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a member of the nucleoporin family. The protein is localized to the nuclear rim and is an essential component of the nuclear pore complex (NPC). All molecules entering or leaving the nucleus either diffuse through or are actively transported by the NPC. Alternate transcriptional splice variants of this gene have been observed but have not been thoroughly characterized.

**Recommended Dilution:** WB,1:100 - 1:500

**Synonyms:** ODG6; NUP84; NPHS11; GAMOS7; Nup107

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 700-800 of human Nup107 (P57740).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.