

CAB20965

Product Information

Product SKU:	CAB20965	Gene ID:	84689	Size:	20uL, 100uL
Clone No:	ARC2954	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	76kDa	Conjugate:	Unconjugated
Calculated MW:	77kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be integral component of membrane.

Recommended Dilution: WB, 1:100 - 1:500

Synonyms: MS4A16; NYD-SP21; MS4A14

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 150-280 of human MS4A14 (Q96JA4).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.