

CAB20967

Product Information

Product SKU:	CAB20967	Gene ID:	54793	Size:	20uL, 100uL
Clone No:	ARC2956	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	45kDa	Conjugate:	Unconjugated
Calculated MW:	43kDa	Isotype:	IgG

Immunogen Information

Background: Enables cullin family protein binding activity; identical protein binding activity; and protein self-association. Predicted to be involved in intracellular signal transduction; protein homooligomerization; and protein ubiquitination.

Recommended Dilution: WB, 1:100 - 1:500

Synonyms: BTBD27; KCTD9

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KCTD9 (Q7L273).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.