

## CAB21088

---

### Product Information

<b>Product SKU:</b>	CAB21088	<b>Gene ID:</b>	140609	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	32kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	35kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** NIMA-related kinases share high amino acid sequence identity with the gene product of the *Aspergillus nidulans* 'never in mitosis A' gene, which controls initiation of mitosis.

**Recommended Dilution:** WB, 1:100 - 1:500

**Synonyms:** -

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 257-301 of human NEK7 (NP\_598001.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.