

CAB21090

Product Information

Product SKU:	CAB21090	Gene ID:	5058/5062/5063	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	61kDa, 66-68kDa	Conjugate:	Unconjugated
Calculated MW:	65kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a family member of serine/threonine p21-activating kinases, known as PAK proteins. These proteins are critical effectors that link RhoGTPases to cytoskeleton reorganization and nuclear signaling, and they serve as targets for the small GTP binding proteins Cdc42 and Rac. This specific family member regulates cell motility and morphology. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2010]

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: PAKalpha; PAK1/PAK2/PAK3

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 350-450 of human PAK1/PAK2/PAK3 (NP_002567.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.