

CAB21095

Product Information

Product SKU:	CAB21095	Gene ID:	1464	Size:	20uL, 100uL
Clone No:	ARC52162	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	250-450kDa	Conjugate:	Unconjugated
Calculated MW:	251kDa	Isotype:	IgG

Immunogen Information

Background: A human melanoma-associated chondroitin sulfate proteoglycan plays a role in stabilizing cell-substratum interactions during early events of melanoma cell spreading on endothelial basement membranes. CSPG4 represents an integral membrane chondroitin sulfate proteoglycan expressed by human malignant melanoma cells.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: NG2; MCSP; MCSPG; MSK16; CSPG4A; HMW-MAA; MEL-CSPG; NG2/CSPG4

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1340-1563 of human NG2/CSPG4 (NP_001888.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH 7.3.