

CAB21116

Product Information

Product SKU:	CAB21116	Gene ID:	4718	Size:	20uL, 100uL
Clone No:	ARC3020	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	14kDa	Conjugate:	Unconjugated
Calculated MW:	14kDa	Isotype:	IgG

Immunogen Information

Background: Involved in mitochondrial respiratory chain complex I assembly. Located in mitochondrion. Part of mitochondrial respiratory chain complex I. Implicated in nuclear type mitochondrial complex I deficiency.

Recommended Dilution: WB,1:500-1:1000

Synonyms: HLC-1; B14.5b; MC1DN36; NADHDH2; CI-B14.5b; NDUFC2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-119 of human NDUFC2 (NP_004540.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.