

## CAB21127

---

### Product Information

<b>Product SKU:</b>	CAB21127	<b>Gene ID:</b>	7812	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC3050	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

---

### Additional Information

<b>Observed MW:</b>	89kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	89kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Enables RNA stem-loop binding activity. Involved in IRES-dependent viral translational initiation; nuclear-transcribed mRNA catabolic process, no-go decay; and stress granule assembly. Located in Golgi apparatus; cytosol; and plasma membrane. Part of CRD-mediated mRNA stability complex.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** UNR; D1S155E; CSDE1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human CSDE1 (NP\_009089.4).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.