

CAB21136

Product Information

Product SKU:	CAB21136	Gene ID:	9456	Size:	20uL, 100uL
Clone No:	ARC3008	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	46kDa	Conjugate:	Unconjugated
Calculated MW:	40kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the homer family of dendritic proteins. Members of this family regulate group 1 metabotropic glutamate receptor function.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: HOMER; SYN47; Ves-1; HOMER1A; HOMER1B; HOMER1C; Homer1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 101-220 of human HOMER1 (NP_004263.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.