

## CAB21143

---

### Product Information

<b>Product SKU:</b>	CAB21143	<b>Gene ID:</b>	11137	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC3017	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	75kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	56kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** The protein encoded by this gene contains several WD-40 repeats and is found mostly in the nucleus. The expression and localization of this protein are cell cycle dependent. Expression of this gene is upregulated in pancreatic adenocarcinoma. Three transcript variants encoding different isoforms have been found for this gene.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** NCLB; IEF-SSP-9502; PWP1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 460-501 of human PWP1 (NP\_008993.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.