

## CAB21144

---

### Product Information

<b>Product SKU:</b>	CAB21144	<b>Gene ID:</b>	11243	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC3037	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	23kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	23kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Enables leucine zipper domain binding activity and transcription coactivator activity. Involved in chromosome segregation. Located in Golgi apparatus; kinetochore; and nucleoplasm. Part of MIS12/MIND type complex. Implicated in bladder carcinoma and urinary bladder cancer.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** PMF1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human PMF1 (NP\_009152.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.