

## CAB21151

---

**Product Information**

<b>Product SKU:</b>	CAB21151	<b>Gene ID:</b>	27297	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC3031	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse

---

**Additional Information**

<b>Observed MW:</b>	20kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	17kDa	<b>Isotype:</b>	IgG

---

**Immunogen Information**

**Background:** This gene encodes a membrane protein that functions as part of a receptor complex for a small neuropeptide that increases intracellular cAMP levels. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.

**Recommended Dilution:** WB, 1:100 - 1:500

**Synonyms:** C17; RCP; RCP9; RPC9; POLR3I; POLR3J; CGRPRCP; CGRP-RCP; CRCP

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human CRCP (NP\_001035737.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.