

CAB21154

Product Information

Product SKU:	CAB21154	Gene ID:	51506	Size:	20uL, 100uL
Clone No:	ARC3001	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	20kDa	Conjugate:	Unconjugated
Calculated MW:	19kDa	Isotype:	IgG

Immunogen Information

Background: UFC1 is an E2-like conjugating enzyme for ubiquitin-fold modifier-1 (UFM1; MIM 610553) (Komatsu et al., 2004 [PubMed 15071506]).

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: NEDSG; HSPC155; UFC1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of Human UFC1 (NP_057490.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.