

CAB21159

Product Information

Product SKU:	CAB21159	Gene ID:	64343	Size:	20uL, 100uL
Clone No:	ARC2987	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	47kDa	Conjugate:	Unconjugated
Calculated MW:	45kDa	Isotype:	IgG

Immunogen Information

Background: AZI2, or NAP1, contributes to the activation of NF κ B (see MIM 164011)-dependent gene expression by activating IKK-related kinases, such as NAK (TBK1; MIM 604834) (Fujita et al., 2003 [PubMed 14560022]).

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: AZ2; NAP1; TILP; AZI2

Purification Method: Affinity purification

Immunogen: Recombinant protein containing a sequence corresponding to amino acids 282-387 of human AZI2 (NP_071906.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.