

## CAB21163

---

### Product Information

<b>Product SKU:</b>	CAB21163	<b>Gene ID:</b>	158158	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC3022	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

---

### Additional Information

<b>Observed MW:</b>	83kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	83kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene is a member of the Rab family of GTPases that are involved in regulation of membrane traffic. The encoded protein contains an N-terminal EF-hand domain, a coiled-coil motif and a C-terminal Rab domain. A potential role as tumor suppressor has been indicated for this gene.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** TSG; RAB45

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human RAB45 (NP\_689786.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.